


# AEC-NASA TECH BRIEF



AEC-NASA Tech Briefs describe innovations resulting from the research and development program of the U.S. AEC or from AEC-NASA interagency efforts. They are issued to encourage commercial application. Tech Briefs are published by NASA and may be purchased, at 15 cents each, from the Clearinghouse for Federal Scientific and Technical Information, Springfield, Virginia 22151.

## Multichannel Analyzers at High Rates of Input



Spectrometer system for use at high rates, employing a commercial multichannel analyzer preceded by a gating system

### The problem:

Improvement of the quality of data from conventional multichannel analyzers that introduce considerable spectral distortion at high rates of input.

### The solution:

When used in conjunction with a well-designed gating system incorporating pole-zero compensation, pileup-rejection, and baseline-restoration, almost any multichannel analyzer having an adequate number of channels can achieve very good resolution at high rates of input. A gating system was added to a conventional Ge(Li) gamma-ray spectrometer. When the input rate was increased from 1,000 to 50,000 pulses per second, the resolution at 1.33 Mev changed only from 3.3 to 3.7 kev, with a spectral shift of 0.08%. These figures should be compared with those for the original ungated spectrometer: resolution increase from 3.3 to 30.0 kev and shift of 3.2%. Furthermore, gating improved the peak-to-continuum ratio by an

order of magnitude at 50,000 pulses per second. At a gross input rate of 100,000 pulses per second obtained from a mixture of  $^{60}\text{Co}$  and  $^{137}\text{Cs}$  sources, the gated system resolved the 1.33-Mev gamma-ray peak to 3.5 kev with a spectral shift of 0.1%.

### How it's done:

The detector is a 13-cm<sup>3</sup>, 17-mm-thick, lithium-drifted planar germanium [Ge(Li)] diode. The input stage of the preamplifier is a cooled FET that is directly coupled to the detector. The coupling networks in the preamplifier are pole-zero compensated, so that the output is virtually a pure, singly differentiated pulse, decaying with a 50-μsec time constant.

The main amplifier, which also incorporates pole-zero compensation, provides single RC integrating and differentiating time constants of 1.2μsec. The amplifier output pulse is passed through a simple baseline-restorer with a time constant selected for

(continued overleaf)

best resolution at medium rates. The signal is delayed 3  $\mu$ sec and presented to the input of a linear gate; it is directly coupled from the restorer through the gating circuit.

The linear gate is opened only when (i) the pulse amplitude falls within the window specified by the single-channel analyzer, (ii) the pileup-rejector fails to detect pulse-shape distortion due to pileup, and (iii) the multichannel analyzer is not busy processing a previous pulse. When the gate opens, the pulse is stretched at its peak amplitude, partially biased off, and amplified to a level suitable for analyzer input. Thus the linear gate serves also as a biased amplifier free of rate and pulse-shape problems. Not only does it restrict the analyzer input to relatively low rates, but also it partially derandomizes the inputs (through the action of the busy signal), thus guaranteeing a substantial quiet period preceding every pulse presented to the analyzer. This quiet time is always greater than 35  $\mu$ sec: sufficient time for the analyzer's adequate recovery.

The usefulness of this spectrometer at high rates is due to the properly designed gating system which allows only a relatively low rate of *selected* pulses to impinge on the input circuitry of the multichannel analyzer. The high-rate features in the gating system and amplifying chain have these principal effects:

1) Pole-zero compensation, particularly in the main amplifier, substantially improves resolution and reduces tailing.

2) Pileup-rejection significantly reduces the rate-contributed continuum at high rates.

3) The baseline-restorer, preceding and directly coupled to the linear gate, is of paramount importance: it all but eliminates spectral shift and produces the major improvement in resolution.

#### Notes:

1. For additional details see W. L. Brown, W. A. Higinbotham, G. L. Miller, R. L. Chase, Eds., *Publ. 1593* (Nat. Acad. Sci.—Nat. Res. Council, Washington, D. C., 1969).

2. Inquiries may be directed to:

Office of Industrial Cooperation  
Argonne National Laboratory  
9700 South Cass Avenue  
Argonne, Illinois 60439  
Reference: B69-10214

Source: S.J. Rudnick and M.G. Strauss  
Electronics Division  
(ARG-10355)

#### Patent status:

Inquiries concerning rights for commercial use of this innovation may be made to:

Mr. George H. Lee, Chief  
Chicago Patent Group  
U.S. Atomic Energy Commission  
Chicago Operations Office  
9800 South Cass Avenue  
Argonne, Illinois 60439